


ACTIVE SURFACE EMG ELECTRODE

Compatible with most EMG systems



Top view

Wide spaced pads with intergrated ground

Wide spaced pads

Narrow spaced pads


Screw-spring adapters for fine wire electrodes

Tool for removing pads

Specifications:

Gain	330
Types of electrode spacing	Narrow - 13/16", Wide - 1 3/8"
Fine wire attachments	Screw-springs
Input impedance	Greater than 100M
CMRR	Typical: 95db
Bandwidth (3db points)	12Hz to 3000Hz

- No skin preparation or gel is required.
- Protected against static electricity.
- Electrode end and connector end are color coordinated for ease of use.
- Quantity discounts available.



B&L Engineering

3002 Dow Ave., Suite 416
Tustin, CA 92780
Web Site: www.bleng.com

Tel: (714) 505-9492
Fax: (714) 505-9493
Email: sales@bleng.com